

ADDENDUM

16 November 2000

THE EFFECT OF CONSECUTIVE BATCH PROCESSING ON PARTICLE SIZE AND THERMAL BEHAVIOR IN A HIGH-SHEAR GRANULATOR

TIMOTHY J. SMITH¹, GARY L. SACKETT¹, TY GOLDSBERRY¹, LIRONG LIU²

¹VECTOR CORPORATION: MARION, IA UNITED STATES ; ²PENWEST PHARMACEUTICALS Co.: PATTERSON, NY UNITED STATES

PURPOSE

To incorporate wet-mass, or high shear, mixer blade power consumption information.

METHODS

In addition to the methods reported in the original poster, the mixer blade power was also monitored. Power trends were overlaid on the "Thermal Trends" graphs previously reported. Revised graphs are shown in the Results section.

Note: Formulations and process parameters are reprinted for reference in Tables 1 and 2.

RESULTS

REVISED CR Formulation A Summary

Parameter	1st Batch	2nd Batch	3rd Batch
D50 (Oven), μm	1182	1056	1054
D50 (FB), μm	1148	1040	1064
ΔT , $^{\circ}\text{C}$	13.3	12.3	12.3
Peak Power, KW	3.6	3.8	3.5
Time to Reach Peak Power, Min	0.3	0.3	0.3

REVISED CR Formulation B Summary

Parameter	1st Batch	2nd Batch	3rd Batch
D50 (Oven), μm	761	794	721
D50 (FB), μm	826	790	767
ΔT , $^{\circ}\text{C}$	13.3	13.0	13.4
Peak Power, KW	2.8	2.6	2.8
Time to Reach Peak Power, Min	2.7	2.3	1.6

REVISED IR Formulation C Summary

Parameter	1st Batch	2nd Batch	3rd Batch
D50 (Oven), μm	1147	1377	1545
D50 (FB), μm	1157	1444	1525
ΔT , $^{\circ}\text{C}$	21.8	22.9	24.0
Peak Power, KW	4.6	4.3	4.6
Time to Reach Peak Power, Min	6.2	6.5	6.0


Table 1 – Formulations

Dry Ingredients	Controlled Release A	Controlled Release B	Immediate Release C
HPMC, K 100 M	30%	--	--
HPMC, K 4 M	--	30%	--
Starch 1500	--	--	15%
MCC, 50M	--	--	30%
Lactose	70%	70%	55%
Total Weight (Kg)	15.70	16.59	18.34
Bulk Density (g/cc)	0.423	0.447	0.489


Table 2 – Processing Parameters

Process Parameters	Controlled Release A	Controlled Release B	Immediate Release C
Pre-Mix Time	3 minutes	3 minutes	3 minutes
Water Infusion Time	10 minutes	8 minutes	8 minutes
Water Added (% , Kg)	35.7% , 8.7 Kg	28.1% , 6.5 Kg	24.7% , 6.0 Kg
Wet Mass Time	5 minutes	7 minutes	7 minutes
Time Between Consecutive Batches	6 to 9 minutes	8 to 9 minutes	6 to 8 minutes


Thermal and Power Trends - CR Formulation A



Thermal and Power Trends - CR Formulation B



Thermal and Power Trends - IR Formulation C



CONCLUSIONS (REGARDING MIXER POWER)

Mixer power trends were slightly different for each of the formulations. Mixer power for CR Formulation A rose almost immediately to peak power during the wet-mass stage and then gradually reduced in magnitude, whereas CR Formulation B rose rapidly, leveled off, and did not reach peak power until 23 - 39% of the wet-mass time was completed. IR Formulation C did not reach peak power until almost the end of the wet-mass time. The small rise and fall of the mixer power at the end of each granulation occurred during the discharge of the product. For the formulations used in this study, temperature rise during wet mass (ΔT) appeared to provide a better endpoint correlation to particle size than peak mixer power.

ACKNOWLEDGMENT

Authors extend their thanks to Mr. Paul Sheskey of the Dow Chemical Company and Colorcon, Global Technical Support for their contribution in this study.